Rema[®]-Sil

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Instructions for use

Mode d'emploi

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DENTAURUM

Rema®-Sil – Technische Daten**Komponente A****Komponente B**

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B 205, T 370, H 300 mm

Rema®-Sil – Technical data

	Component A	Component B
Mixing ratio:	1	1
Viscosity: at 23 °C/73.4 °F mPas:	1000	1000
Density kg/m³:	1060	1070
Working time:	4 – 5 minutes	
Setting time:	20 – 30 minutes	

These instructions are based on the results of meticulous tests. We pass them on without further liability. Further individual work-safety regulations should be complied with strictly.

Vulcanisate:

Mechanical properties after vulcanization at 23 °C/73.4 °F

Measured with a 6 mm thick test sample:

Shore A hardness (ASTM D 2240) ca. 9-10 units

Measured with a 2 mm thick test sample:

Tensile strength (AFNOR NF T 46 002) 1,2 – 1,5 MPa

Extension at rupture (AFNOR NF T 46 002) 450 – 550 %

Temperature stability: ca. –70 to + 150 °C/ –158 to + 302 °F

Inhibition

If Rema®-Sil comes into contact with the following substances, vulcanization may be delayed or prevented altogether:

- condensation cross-linked silicones
- PVC stabilizers
- sulferous natural rubber
- to a lesser extent, modeling clay
- Oils (milling oil)

Health and Safety Hazards

Avoid swallowing as well as prolonged or repeated skin contact with Rema®-Sil components A and B.

In case of accidental skin contact, wash thoroughly with soap and water.

In case of contact with the eyes wash the eyes thoroughly for 15 minutes.

Contact an opthamologist immediately.

Rema®-Sil – Instructions for use


is an additional cross-linked two component silicone for the production of duplicate models in dental prosthetics. Mixing ratio 1 : 1. Thin. Translucent blue. Easy to mix with any vacuum mixing unit.

Uses

Suitable for making moulds that may be poured with: Plaster such as Rapidur, Quartzous investment materials such as Rema® Star, Rema® Exakt, Castorit®-super, Castorit®-super C, Methyl methacrylate based acrylics such as Orthocryl® Rapid.

Important! Storage and Working Temperatures.

Store at room temperature of 20 °C (68 °F) to 23 °C (73 °F). If the Rema®-Sil has been stored at a lower temperature it must be brought to room temperature before usage. Due to its lower thermal conductivity, 1 to 2 days storage time at a higher temperature will be required. The time needed to raise the temperature will depend upon the storage temperature. Disregarding these temperature requirements will result in delayed or incomplete setting of the material. During cool seasons wait an additional 20 minutes before removing the master cast. This will avoid the risk of mold deformation due to the silicone not being completely set.

The setting time of silicone depends on the temperature.

Higher than recommended storage or room temperature will accelerate the setting of the silicone.

Note: In order to produce perfect models

Pouring models with other products may result in imperfect surface characteristics due to product incompatibility. We therefore recommend preliminary testing before proceeding with actual work.

Duplication**The mastermodel**

Use dry, clean (steam cleaning) plaster models. Block-out undercut areas with Gumex. Place a small roll of Gumex on the distal end of the model which, when removed will build a cavity. This can be used to deflask the model with compressed air. Spray the gypsum surface with Septisol release agent and blow dry immediately.

Dispensing Rema®-Sil from the 5 kg Component A+B Containers**a) Drain cock**

Lay the container horizontally on a stand or shelf with the drain cock facing downwards. In order to assure adequate silicone flow bore a small, 5 mm, hole to the right or left of the grip.

b) In the Silimix unit

When the button is pressed, the silicone components are forced in a ratio of 1 : 1 through a nozzle, mixed in the nozzle and introduced directly into the flask. When the button is released, the unit switches off. To guarantee good subsequent delivery of the silicone, a hole of approx. Ø 5 mm is drilled at the highest point in each component's container.

Using Rema®-Sil

1. Automatic mixing of parts A and B in Silimix unit.
2. Weigh out 1 parts of component A and 1 of component B of Rema®-Sil.

We recommend the Soehnle scale Type 800 900 with electronic zero-setting system (1 g).

Place a clean mixing beaker from the vacuum mixing unit and set the display to "0".

Dispense the desired amount of Rema®-Sil, eg. 50 g = 1 part.

Reset the display to "0" and dispense component B (blue bottle), eg. 50 g = 1 part.

Important! Pour slowly since the read-out lags behind the actual weight.

Mix the two components in a mixing unit for 60 seconds with vacuum. Pour the mixture into the flask from a height of about 30 cm (between the flask and the mixing beaker).

Any rests remaining in the beaker can be removed after 20 – 30 minutes.

(Removal can be made faster by heating the mixing beaker).

No pressure is necessary in order to cure the silicone. When mixed properly in a vacuum, and if the model has been prepared correctly, no bubbles will form on the surface of the mould. If several consecutive operations are to be carried out, the use of two beakers is recommended. Setting time: 40 minutes. Removal of the model before the fully elapsed setting time may result in deformations in the silicone since the proper Shore hardness has not been attained.

Removing the Model

Remove the small role of Gumex from the model with a small instrument. Carefully blow out the model with compressed air.

Preparing the Mould

In order to assure a smooth flow of the filling material spray the mould with Lubrofilm® degreasing agent.

Blow dry immediately. Important! A moist silicone surface will result in an imperfect model surface.

Important!

If the silicone surface is not treated with degreaser, the filling material will not flow smoothly and bubbles may result.

Pouring the Mould

Investment materials with high tensile strength such as Rema® Star are recommended. Rema® Star does not have to be hardened. Dry out the models in a hot air drying cabinet for 20 minutes (depending upon the size of the model) at 70 °C (158 °F). **Undried** models should be waxed-up immediately, invested, as well as burned-out and cast on the same day.

Re-use of Hardened Silicone

Rema®-Sil is irreversible, however, chopped-up pieces of silicone (the smaller the better) can be used as a filler in combination with newly mixed material. There is absolutely no effect on the accuracy of fit.

Disposal

Rema®-Sil can be disposed of when set, with domestic rubbish.

Shelf life

12 months provided it is stored under dry conditions.

Trouble shooting guide

Problems	Probable causes	Solutions
1. Rough areas on model	The two components A & B not mixed correctly.	Mixing for 60 seconds until a uniform blue color is achieved.
2. Tracky areas on model. Silicone not set.	Grease and insulating agent residues left on model. A & B components not mixed correctly. Lubricating oil on the gypsum surface. Unsuitable isolating agent. Residue from light curing materials. Unsuitable gypsum hardener.	Clean models with mild detergent in water, then with steam jet; spray with Septisol and blow dry. When using mixing units, dispense a little silicone separately from the nozzle, and only then fill the flask. 25–30 minutes hardening time in pressure pot. Clean the model with a diluted mild detergent and a sponge, steam clean thoroughly. Construct the milling model. Clean the model with a diluted mild detergent and a sponge and steam clean thoroughly. Remove the residue with alcohol and use an alternative material, e. g. waxes. Use an alternative material. Do not harden work to be duplicated.
3. Pock-marks and flaking investment on the model.	Model released from the silicone too early.	40 minutes hardening time for the investment material.
4. Soft models.	Heavily sprayed with Lubrofilm® and damp mould filled with investment material.	Spray Lubrofilm® into the silicone out mould and blow dry immediately.
5. Warped castings.	Incorrect mixing ratio. Master model does not release smoothly from silicone. Silicone loosened from flask. Difficulty removing master Models. Button taped on removal from investment.	Mix A & B 1 : 1. Block-out undercut areas thoroughly. Overfill opening of flask for better retention. Apply a little Gumex, to model socket for better release.
6. Rema®-Sil does not set or only partly sets.	Oil substances, e. g. milling oil, silicones contaminated by condensation, natural rubber containing sulphur, stabilizers and PVC and to a lesser extent standard plastics were used.	Avoid these materials listed.

Recommended Duplicating Systems

Best results are achieved with our Rema®-Star model casting systems. See special brochure. Dentaurem's Siliform transparent duplicating flasks in 3 sizes can also be used.

Delivery Program

Rema®-Sil

		Ord.-Nr.
Component A	1 kg	108-700
Component B	1 kg	108-701
Component A	5 kg	108-710
Component B	5 kg	108-711

Septisol, die-lubricating agent for plaster/silicone and plaster/plaster	
100 ml	108-720
1000 ml	108-721

Lubrofilm®, wax and silicone surface tension reducing agent.	
100 ml	112-050
1000 ml	112-051

Gumex, lining compound, sticking to the model, 10 sticks of 13 cm, Ø 7 mm	168-015
---	----------------

Drain cock for 5 kg container	108-716
Measuring and mixing beaker	108-715

Neo-Star Multi-Function Flask

Special duplicating flask featuring: Savings in silicone and investment materials, duplicating, repositioning, investing, model-base mounting, articulation, magnetic fixing.

Standard Neo-Star flask/red	165-800
Large Neo-Star flask/green	165-810

Siliform Duplicating Flask

3 sizes, transparent, plastic.

Siliform set with 3 flasks **127-318**
small, medium, large

1 Top section, small **127-315**

1 Top section, medium **127-316**

1 Top section, large **127-317**

1 Base for all sizes **127-309**

Rema® Star Investment Material

super hard, silk-smooth model surface.
No hardening necessary.

80 x 250 g Packs = 20 kg **105-700**

112 x 180 g Packs = 20 kg **105-705**

1000 ml mixing fluid **105-501**
(3.2 liters are required to mix 20 kg).

Airvac – Vacuum Mixing Unit

for silicones, investment materials and plaster.
Unit includes 2 stirring beakers

Wall unit 230 volts **095-070**

L 340, W 350, H 285 mm

Bench unit 230 volts **095-060**

L 340, W 350, H 555 mm

Silimix

Silicone measuring and mixing unit

For two-component silicones with a mixing
ratio of 1 : 1

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